

Elastic Ankle Support

Product description

Designed to provide light support following stressed related ankle injuries. It does not restrict movements and may be used during physical activity.

Characteristics

- Even and comfortable compression
- Light support
- Easy adjustment
- Washable and reusable
- Made of 77% nylon, 23% elastic
- Colour: black

Made of: 77% nylon, 23% elastic

Indications

- Swelling
- Ankle instability

Available sizes

| CODE | ANKLE CIRCUMFERENCE | SIZE |
|------|--------------------------|--------|
| 9282 | 17 - 21 cm (6.5" - 8") | SMALL |
| 9283 | 21 - 24 cm (8" - 9.5") | MEDIUM |
| 9284 | 24 - 29 cm (9.5" - 11") | LARGE |
| 9285 | 29 - 32 cm (11" - 12.5") | XLARGE |

To determine proper size, measure the circumference of the ankle. If the result is between two sizes, choose the largest.

Adjustment

Pull-up the ankle sleeve and adjust accordingly.

Cleaning

Wash by hand with mild soap and air-dry.

Caution: If pain should persist, or if you wish to know more about how to use this product, consult a health specialist.

