

## Elastic Knee Support

### Product description

Designed to provide light support following minor knee injuries. The Elastic Knee Support does not restrict movements and may be used during physical activity.

### Characteristics

- Offers a light support after a trauma to the knee
- Tension reduced at patella level
- Light support
- Easy adjustment
- Made of 77% nylon and 23% elastic
- Colour: black

**Made of:** 77% nylon, 23% elastic

### Indications

- to wear during activities which lead to discomfort
- during physical activity to reduce the risk of injuries

### Available sizes

CODE	KNEE CIRCUMFERENCE	SIZE
9316	35 - 38 cm (14" - 15")	SMALL
9317	38 - 43 cm (15" - 17")	MEDIUM
9318	43 - 51 cm (17" - 20")	LARGE
9319	51 - 58 cm (20" - 23")	XLARGE

To determine proper size, measure the circumference at the mid-section of the knee joint. If the result is between two sizes, choose the largest.

### Adjustment

Pull-up the knee sleeve and adjust accordingly.

### Cleaning

Wash by hand with mild soap and air-dry.

**Caution:** If pain should persist, or if you wish to know more about how to use this product, consult a health specialist.

