

Arm cast support

Product description

Designed to reduce cervical fatigue by transferring the weight of injured hand over the opposite shoulder. May be used for right or left arm. Efficient for the entire convalescence period of the patient.

Characteristics

- Solid polypropylene strap for easy adjustment
- Padding on neck line can be adjusted for optimal comfort
- Made of polyurethane foam laminated with brushed nylon; terry cloth interior
- Colour child: grey with white piping
- Colour adult: grey with black piping

Made of: Laminated brushed nylon on the exterior on synthetic foam and terry cloth interior.

Indications

- Injury to the arm, wrist or hand

Available sizes

CODE	MODEL
7679	KIDS
7680	ADULT

Adjustment

Place the support under the injured arm. Adjust the felt pad on the neck line on the opposite shoulder for an optimal comfort.

Cleaning

Wash by hand with mild soap and air-dry.

