

FORMEDICA™ Leg and Knee Immobilizer

Product description

Recommended for sprains, dislocations, minor fractures, various arthropathies and before or after surgery.

All immobilizers packages include a hygienic hosing.

Indications

Recommended in case of sprains, dislocations, simple fractures, and before or after surgery.

Characteristics

- Five removable and malleable aluminum stays: one on each side and three in the back
- Made of polyurethane foam laminated with cotton; terry cloth interior
- Velcro®-type closures

Available sizes

HEIGHT 49cm (19")	HEIGHT 56cm (22")	HEIGHT 64cm (25")	COLOUR	ADJUSTMENT	SIZE
7376	8401	8409	GREY	38 – 43cm (15" – 17")	SMALL
7377	8402	8410	GREEN	43 – 49cm (17" – 19")	MEDIUM
7378	8403	8411	BLUE	49 – 54cm (19" – 21")	LARGE
7379	8404	8412	GREY	54 – 59cm (21" – 23")	XLARGE



Cleaning

Remove stays. Wash on delicate cycle and air-dry.