

Knee Support

Product description

Designed to help reduce the risk of injuries during or after sports activities. The Knee Support provides compression and heat to support weak knee muscles and joints following strains, sprains and other stress-related injuries. Neoprene retains body heat, increases blood circulation and helps keep muscles supple. Remains effective when wet and does not absorb water.

Characteristics

- Sleeve support
- Remains effective when wet (does not absorb water)
- Provides therapeutic heat and compression
- Neoprene covered with non-brushed nylon
- Colour: black

Made of: Neoprene covered with non-brushed nylon.

Indications

- Strains
- Sprain
- Post-traumatic usage
- Sports

Available sizes

CODE	KNEE CIRCUMFERENCE	SIZE
9437	33 - 36 cm (13" - 14")	SMALL
9438	36 - 38 cm (14" - 15")	MEDIUM
9439	38 - 43 cm (15" - 17")	LARGE
9440	43 - 51 cm (17" - 20")	XLARGE

To determine proper size, measure the circumference at the mid-section of the knee joint. If the result is between two sizes, choose the largest.

Adjustment

Pull up the knee sleeve and make sure the opening on the knee support is around the patella.

Cleaning

Wash by hand with mild soap and air-dry.

Caution: If pain should persist, or if you wish to know more about how to use this product, consult a health specialist.

