

Aluminum Hand and Wrist Splint

Product description

Formedica's aluminum splints allow to immobilize temporarily fingers and wrist in case of sprain or fracture. Their aluminum side provides the necessary rigidity and their foam side ensures patient's comfort.



Indications

- To be used in cases of wrist or radius fractures.
- Injuries in various sports (basketball, volleyball, etc ...)

Characteristics

- Easily malleable
- Radiolucent
- Light aluminum covered with polyurethane foam
- Non-allergenic
- Can be easily cut to desired length

Available sizes

CODE DOZEN	DIMENSIONS
127926	5 cm x 46 cm (2" x 18")
127927	7,5 cm x 41 cm (3" x 16")